

Chemical Polarimeter Cells

(Order Code CELLS-POL)

Chemical Polarimeter Cells are a pack of four cells for use with the Vernier Chemical Polarimeter. They are the same sample cells as the one that ships with the Chemical Polarimeter. Each cell is a flat-bottom glass tube with a 20 mm inner diameter and is 150 mm tall. A ruler printed on the tube displays cm markings up to 12 cm high to measure sample height. They are designed to fit snugly in the Chemical Polarimeter to reduce inconsistencies between runs. There is also an easy-pour spout for transferring solutions safely and efficiently.



Vernier Software & Technology

13979 S.W. Millikan Way • Beaverton, OR 97005-2886

Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440

info@vernier.com • www.vernier.com