

Quartz Cuvettes

(Order Code CUV-QUARTZ)

The Quartz Cuvettes are a pack of two Spectrosil® Quartz cuvettes for use with UV-VIS spectrophotometers. They are the same pack of cuvettes that ship with the Vernier UV-VIS Spectrophotometer. Each cuvette is a standard rectangular cuvette with two polished sides.



- Spectral range: 170–2700 nm
- Path length: 10 mm
- Exterior dimensions: 12.5 mm (width), 12.5 mm (length), and 45 mm (height)
- Interior dimensions: 10 mm (width) and 10 mm (length)
- Nominal volume: 3.5 mL

Vernier Software & Technology

13979 S.W. Millikan Way • Beaverton, OR 97005-2886

Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440

info@vernier.com • www.vernier.com